

Ack1 (Phospho-Tyr857+Tyr858) Antibody

Catalog No: #12658



Package Size: #12658-1 50ul #12658-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Ack1 (Phospho-Tyr857+Tyr858) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu
Specificity	Ack1 (Phospho-Tyr857+Tyr858) Antibody detects endogenous levels of Ack1 only when phosphorylated at Tyr857+Tyr858
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Ack1 (Phospho-Tyr857+Tyr858)
Target Name	Ack1
Modification	Phospho
Other Names	TNK2, ACK, Activated p21cdc42Hs kinase, Ack kinase, ACK-1, ACK1, Activated CDC42 kinase 1, p21cdc42Hs
Accession No.	Swiss-Prot#: Q07912/Q9UJM3NCBI Gene ID: 10188/54206
Target Species	human
Calculated MW	115kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of SH-SY5Y, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.