

Vasopressin V1b receptor

Catalog No: #AP89265



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Vasopressin V1b receptor |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | V1bR,AVPR V1b,AVPR V3,Antidiuretic hormone receptor 1b,Vasopressin V3 receptor,V1br |
| Accession No. | Uniprot ID:Q9WU02 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.