

Vascular endothelial growth factor receptor 3

Catalog No: #AP89298



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor receptor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	VEGFR-3, Fms-like tyrosine kinase 4, FLT-4, Tyrosine-protein kinase receptor FLT4, Flt-4, Vegfr3
Accession No.	Uniprot ID:Q91ZT1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.