

ABI1 (Phospho-Tyr212) Antibody

Catalog No: #12670

Package Size: #12670-1 50ul #12670-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABI1 (Phospho-Tyr212) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ABI1 (Phospho-Tyr212) Antibody detects endogenous levels of ABI1 only when phosphorylated at Tyr212
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ABI1 (Phospho-Tyr212)
Target Name	ABI1
Modification	Phospho
Other Names	ABI1, Abelson interactor 1, Abl-interactor protein 1 long, ABI-1, Abl-interactor 1, Abl interactor 1, Abl-binding protein 4, ABLBP4, NAP1BP, Nap1 binding protein, SSH3BP1, Nap1-binding protein, E3B1, Eps8-binding protein, Interactor protein AbIBP4, SSH3BP
Accession No.	Swiss-Prot#: Q8IZP0NCBI Gene ID: 10006
Target Species	human
Calculated MW	55kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.