

## Chicken ALCAM/CD166 antigen ELISA Kit

Catalog No: #EK19156

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken ALCAM/CD166 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken CD166 antigen
Target Name	CD166 antigen
Other Names	Activated leukocyte cell adhesion molecule, BEN glycoprotein, Protein DM-GRASP, Protein JC7, SC1 glycoprotein, CD166
Accession No.	P42292

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.15ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ALCAM

Note: This product is for in vitro research use only and is not intended for use in humans or animals.