

Chicken ACE/Angiotensin-converting enzyme ELISA Kit

Catalog No: #EK19159



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chicken ACE/Angiotensin-converting enzyme ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Angiotensin-converting enzyme
Target Name	Angiotensin-converting enzyme
Other Names	ACE, Dipeptidyl carboxypeptidase I, Kininase II, DCP1
Accession No.	Q10751

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.15ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ACE

Note: This product is for in vitro research use only and is not intended for use in humans or animals.