

Pig ANG/Angiogenin ELISA Kit

Catalog No: #EK19162



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pig ANG/Angiogenin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Angiogenin
Target Name	Angiogenin
Other Names	Ribonuclease 5, RNase 5
Accession No.	P31346

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 21pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANG

Note: This product is for in vitro research use only and is not intended for use in humans or animals.