

Bovine ANGPT1/Angiopietin-1 ELISA Kit

Catalog No: #EK19163



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine ANGPT1/Angiopietin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Angiopietin-1
Target Name	Angiopietin-1
Other Names	ANG-1,ANG1
Accession No.	O18920

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.15ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ANGPT1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.