

Bovine BDNF/Brain-derived neurotrophic factor ELISA Kit

Catalog No: #EK19169



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine BDNF/Brain-derived neurotrophic factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Brain-derived neurotrophic factor
Target Name	Brain-derived neurotrophic factor
Other Names	BDNF
Accession No.	Q95106

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BDNF

Note: This product is for in vitro research use only and is not intended for use in humans or animals.