

# Chicken BMP2/Bone morphogenetic protein 2 ELISA Kit

Catalog No: #EK19172



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken BMP2/Bone morphogenetic protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Bone morphogenetic protein 2
Target Name	Bone morphogenetic protein 2
Other Names	BMP-2,BMP-2
Accession No.	Q90751

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:BMP2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.