

# Bovine BMP4/Bone morphogenetic protein 4 ELISA Kit

Catalog No: #EK19174



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine BMP4/Bone morphogenetic protein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Bone morphogenetic protein 4
Target Name	Bone morphogenetic protein 4
Other Names	BMP-4
Accession No.	Q2KJH1

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 37pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: BMP4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.