

Chicken BMP4/Bone morphogenetic protein 4 ELISA Kit

Catalog No: #EK19175



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chicken BMP4/Bone morphogenetic protein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Bone morphogenetic protein 4
Target Name	Bone morphogenetic protein 4
Other Names	BMP-4,BMP-4
Accession No.	Q90752

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 5.4pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BMP4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.