

Rat Bmp4/Bone morphogenetic protein 4 ELISA Kit

Catalog No: #EK19176



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Bmp4/Bone morphogenetic protein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bone morphogenetic protein 4
Target Name	Bone morphogenetic protein 4
Other Names	BMP-4,Bone morphogenetic protein 2B,BMP-2B,Bmp-4,Dvr-4
Accession No.	Q06826

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bmp4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.