

# Human BMPR2/Bone morphogenetic protein receptor type-2 ELISA Kit

Catalog No: #EK19178



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human BMPR2/Bone morphogenetic protein receptor type-2 ELISA Kit                          |
| Brief Description  | ELISA Kits  |
| Applications       | ELISA   |
| Species Reactivity | Hu  |
| Specificity        | Natural and recombinant Human Bone morphogenetic protein receptor type-2                  |
| Target Name        | Bone morphogenetic protein receptor type-2  |
| Other Names        | BMP type-2 receptor,Bone morphogenetic protein receptor type II,BMP type II receptor,PPH1 |
| Accession No.      | Q13873  |

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:BMPR2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.