

Mouse Cdh1/Cadherin-1 ELISA Kit

Catalog No: #EK19183



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cdh1/Cadherin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cadherin-1
Target Name	Cadherin-1
Other Names	ARC-1, Epithelial cadherin, E-cadherin, Uvomorulin, CD324
Accession No.	P09803

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.1pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdh1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.