

Rat Cdh1/Cadherin-1 ELISA Kit

Catalog No: #EK19184



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cdh1/Cadherin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cadherin-1
Target Name	Cadherin-1
Other Names	Epithelial cadherin,E-cadherin,Uvomorulin,CD324
Accession No.	Q9R0T4

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.084ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cdh1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.