

# Chicken CNTF/Ciliary neurotrophic factor ELISA Kit

Catalog No: #EK19185



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chicken CNTF/Ciliary neurotrophic factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Chicken
Specificity	Natural and recombinant Chicken Ciliary neurotrophic factor
Target Name	Ciliary neurotrophic factor
Other Names	CNTF, Growth-promoting activity, GPA, CNF, GPA
Accession No.	Q02011

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CNTF

Note: This product is for in vitro research use only and is not intended for use in humans or animals.