

Bovine CORO1A/Coronin-1A ELISA Kit

Catalog No: #EK19188



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine CORO1A/Coronin-1A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Coronin-1A
Target Name	Coronin-1A
Other Names	Coronin-like protein A,Clipin-A,Coronin-like protein p57,Tryptophan aspartate-containing coat protein,TACO,CORO1
Accession No.	Q92176

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CORO1A

Note: This product is for in vitro research use only and is not intended for use in humans or animals.