

# Bovine CSF3/Granulocyte colony-stimulating factor ELISA Kit

Catalog No: #EK19191



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CSF3/Granulocyte colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Granulocyte colony-stimulating factor
Target Name	Granulocyte colony-stimulating factor
Other Names	G-CSF,GCSF
Accession No.	P35833

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CSF3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.