

Mouse Csf3/Granulocyte colony-stimulating factor ELISA Kit

Catalog No: #EK19193



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Csf3/Granulocyte colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Granulocyte colony-stimulating factor
Target Name	Granulocyte colony-stimulating factor
Other Names	G-CSF, Csf3
Accession No.	P09920

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.64pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Csf3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.