

Pig CSF3/Granulocyte colony-stimulating factor ELISA Kit

Catalog No: #EK19194



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pig CSF3/Granulocyte colony-stimulating factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Granulocyte colony-stimulating factor
Target Name	Granulocyte colony-stimulating factor
Other Names	G-CSF
Accession No.	O02837

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CSF3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.