

# Pig CCL8/C-C motif chemokine 8 ELISA Kit

Catalog No: #EK19201



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig CCL8/C-C motif chemokine 8 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig C-C motif chemokine 8
Target Name	C-C motif chemokine 8
Other Names	Monocyte chemoattractant protein 2, Monocyte chemotactic protein 2, MCP-2, Small-inducible cytokine A8, MCP2, SCYA8
Accession No.	P49873

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CCL8

Note: This product is for in vitro research use only and is not intended for use in humans or animals.