

Mouse Csf1/Macrophage colony-stimulating factor 1 ELISA Kit

Catalog No: #EK19204



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Csf1/Macrophage colony-stimulating factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Macrophage colony-stimulating factor 1
Target Name	Macrophage colony-stimulating factor 1
Other Names	CSF-1,Csfm
Accession No.	P07141

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.09ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Csf1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.