

Mouse Ccl22/C-C motif chemokine 22 ELISA Kit

Catalog No: #EK19206



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ccl22/C-C motif chemokine 22 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C motif chemokine 22
Target Name	C-C motif chemokine 22
Other Names	Activated B and dendritic cell-derived,CC chemokine ABCD-1,Small-inducible cytokine A22,Abcd1,Scya22
Accession No.	O88430

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.87pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ccl22

Note: This product is for in vitro research use only and is not intended for use in humans or animals.