

# Bovine CCL3/C-C motif chemokine 3 ELISA Kit

Catalog No: #EK19207



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CCL3/C-C motif chemokine 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine C-C motif chemokine 3
Target Name	C-C motif chemokine 3
Other Names	Macrophage inflammatory protein 1-alpha, MIP-1-alpha, Small-inducible cytokine A3, SCYA3
Accession No.	Q8SQA6

## Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CCL3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.