

Chicken CCL4/C-C motif chemokine 4 homolog ELISA Kit

Catalog No: #EK19210



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Chicken CCL4/C-C motif chemokine 4 homolog ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Chicken |
| Specificity | Natural and recombinant Chicken C-C motif chemokine 4 homolog |
| Target Name | C-C motif chemokine 4 homolog |
| Other Names | Macrophage inflammatory protein 1-beta homolog,MIP-1-beta,Small-inducible cytokine A4 homolog,SCYA4 |
| Accession No. | Q90826 |

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CCL4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.