

Bovine CCL20/C-C motif chemokine 20 ELISA Kit

Catalog No: #EK19214



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine CCL20/C-C motif chemokine 20 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine C-C motif chemokine 20
Target Name	C-C motif chemokine 20
Other Names	Macrophage inflammatory protein 3 alpha, MIP-3-alpha, Small-inducible cytokine A20, MIP3A, SCYA20
Accession No.	Q8SQB1

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL20

Note: This product is for in vitro research use only and is not intended for use in humans or animals.