

Rat Ccl20/C-C motif chemokine 20 ELISA Kit

Catalog No: #EK19216

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ccl20/C-C motif chemokine 20 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat C-C motif chemokine 20
Target Name	C-C motif chemokine 20
Other Names	Beta-chemokine exodus-1,CC chemokine LARC,CC chemokine ST38,Liver and activation-regulated chemokine,Macrophage inflammatory protein 3 alpha,MIP-3-alpha,Small-inducible cytokine A20,Scya20,St38
Accession No.	P97884

Application Details

Assay Type: Sandwich Detect Range: 15.6-1000pg/mL Sensitivity: 7.92pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ccl20

Note: This product is for in vitro research use only and is not intended for use in humans or animals.