

# Dog CCL5/C-C motif chemokine 5 ELISA Kit

Catalog No: #EK19219



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog CCL5/C-C motif chemokine 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog C-C motif chemokine 5
Target Name	C-C motif chemokine 5
Other Names	Small-inducible cytokine A5,T-cell-specific protein RANTES,SCYA5
Accession No.	Q8HYS0

## Application Details

Assay Type: Sandwich Detect Range:39-2500pg/mL Sensitivity:19.5pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CCL5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.