

Pig CCL5/C-C motif chemokine 5 ELISA Kit

Catalog No: #EK19221



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pig CCL5/C-C motif chemokine 5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig C-C motif chemokine 5
Target Name	C-C motif chemokine 5
Other Names	SIS-delta, Small-inducible cytokine A5, T-cell-specific protein RANTES, SCYA5
Accession No.	Q29288

Application Details

Assay Type: Sandwich Detect Range: 39-2500 pg/mL Sensitivity: 19.5 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.