

Human COL28A1/Collagen alpha-1(XXVIII) chain ELISA Kit

Catalog No: #EK19222



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human COL28A1/Collagen alpha-1(XXVIII) chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Collagen alpha-1(XXVIII) chain
Target Name	Collagen alpha-1(XXVIII) chain
Other Names	COL28
Accession No.	Q2UY09

Application Details

Assay Type: Sandwich Detect Range: 0.156-10nmol/L Sensitivity: 0.032nmol/L Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: COL28A1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.