

# Dog CD40LG/CD40 ligand ELISA Kit

Catalog No: #EK19225



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dog CD40LG/CD40 ligand ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog CD40 ligand
Target Name	CD40 ligand
Other Names	CD40-L, Tumor necrosis factor ligand superfamily member 5, CD40L, TNFSF5, CD154
Accession No.	O97626

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.15ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CD40LG

Note: This product is for in vitro research use only and is not intended for use in humans or animals.