

# Bovine BDH1/D-beta-hydroxybutyrate dehydrogenase, mitochondrial ELISA Kit

Catalog No: #EK19231



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine BDH1/D-beta-hydroxybutyrate dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine D-beta-hydroxybutyrate dehydrogenase, mitochondrial
Target Name	D-beta-hydroxybutyrate dehydrogenase, mitochondrial
Other Names	3-hydroxybutyrate dehydrogenase,BDH,BDH
Accession No.	Q02337

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:BDH1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.