

## Human CRYAB/Alpha-crystallin B chain ELISA Kit

Catalog No: #EK19233

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CRYAB/Alpha-crystallin B chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alpha-crystallin B chain
Target Name	Alpha-crystallin B chain
Other Names	Alpha(B)-crystallin, Heat shock protein beta-5, HspB5, Renal carcinoma antigen NY-REN-27, Rosenthal fiber component, CRYA2, HSPB5
Accession No.	P02511

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.089ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CRYAB

Note: This product is for in vitro research use only and is not intended for use in humans or animals.