

Mouse Cryab/Alpha-crystallin B chain ELISA Kit

Catalog No: #EK19234



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cryab/Alpha-crystallin B chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Alpha-crystallin B chain
Target Name	Alpha-crystallin B chain
Other Names	Alpha(B)-crystallin,P23,Crya2
Accession No.	P23927

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/ml Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cryab

Note: This product is for in vitro research use only and is not intended for use in humans or animals.