

# Rat Cryab/Alpha-crystallin B chain ELISA Kit

Catalog No: #EK19235

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cryab/Alpha-crystallin B chain ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Alpha-crystallin B chain
Target Name	Alpha-crystallin B chain
Other Names	Alpha(B)-crystallin
Accession No.	P23928

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10 ng/mL Sensitivity: 0.078 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100 µl Assay Time: 3-5 h Detection wavelength: 450 nm gene\_full\_name: Cryab

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.