

Mouse Cav1/Caveolin-1 ELISA Kit

Catalog No: #EK19236



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cav1/Caveolin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Caveolin-1
Target Name	Caveolin-1
Other Names	Cav
Accession No.	P49817

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cav1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.