

Rat Adm/ADM ELISA Kit

Catalog No: #EK19239



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Rat Adm/ADM ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat ADM
Target Name	ADM
Accession No.	P43145

Application Details

Assay Type: SandwichDetect Range:7.8-500pg/mL Sensitivity:4.17pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Adm

Note: This product is for in vitro research use only and is not intended for use in humans or animals.