

Pig AMH/Muellerian-inhibiting factor ELISA Kit

Catalog No: #EK19240



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pig AMH/Muellerian-inhibiting factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Muellerian-inhibiting factor
Target Name	Muellerian-inhibiting factor
Other Names	Anti-Muellerian hormone,AMH,Muellerian-inhibiting substance,MIS
Accession No.	P79295

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.091ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:AMH

Note: This product is for in vitro research use only and is not intended for use in humans or animals.