

# Human CDO1/Cysteine dioxygenase type 1 ELISA Kit

Catalog No: #EK19241



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CDO1/Cysteine dioxygenase type 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cysteine dioxygenase type 1
Target Name	Cysteine dioxygenase type 1
Other Names	Cysteine dioxygenase type I,CDO
Accession No.	Q16878

## Application Details

Assay Type: SandwichDetect Range:0.5-32ng/mL Sensitivity:0.23ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:CDO1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.