

Rat Alpi/Intestinal-type alkaline phosphatase 1 ELISA Kit

Catalog No: #EK19244



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Alpi/Intestinal-type alkaline phosphatase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Intestinal-type alkaline phosphatase 1
Target Name	Intestinal-type alkaline phosphatase 1
Other Names	IAP-1, Intestinal alkaline phosphatase I, IAP-I
Accession No.	P15693

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.88ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Alpi

Note: This product is for in vitro research use only and is not intended for use in humans or animals.