

Mouse Crhr2/Corticotropin-releasing factor receptor 2 ELISA Kit

Catalog No: #EK19245

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Crhr2/Corticotropin-releasing factor receptor 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Corticotropin-releasing factor receptor 2
Target Name	Corticotropin-releasing factor receptor 2
Other Names	CRF-R-2, CRF-RB, Corticotropin-releasing hormone receptor 2, CRH-R-2, Crf2r
Accession No.	Q60748

Application Details

Assay Type: Sandwich Detect Range: 0.312-20 ng/mL Sensitivity: 0.161 ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100 µL Assay Time: 3-5 h Detection wavelength: 450 nm gene_full_name: Crhr2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.