

Mouse Acat1/Acetyl-CoA acetyltransferase, mitochondrial ELISA Kit

Catalog No: #EK19249



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Acat1/Acetyl-CoA acetyltransferase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Acetyl-CoA acetyltransferase, mitochondrial
Target Name	Acetyl-CoA acetyltransferase, mitochondrial
Other Names	Acetoacetyl-CoA thiolase
Accession No.	Q8QZT1

Application Details

Assay Type: Sandwich Detect Range: 7.8-500pg/mL Sensitivity: 4.18pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Acat1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.