

# Human ATOH1/Protein atonal homolog 1 ELISA Kit

Catalog No: #EK19252



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ATOH1/Protein atonal homolog 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Protein atonal homolog 1
Target Name	Protein atonal homolog 1
Other Names	Class A basic helix-loop-helix protein 14,bHLHa14,Helix-loop-helix protein hATH-1,hATH1,ATH1,BHLHA14
Accession No.	Q92858

## Application Details

Assay Type: SandwichDetect Range:0.78-50ng/mL Sensitivity:0.33ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:ATOH1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.