

# Mouse Arsa/Arylsulfatase A ELISA Kit

Catalog No: #EK19254



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Mouse Arsa/Arylsulfatase A ELISA Kit          |
| Brief Description  | ELISA Kits                                    |
| Applications       | ELISA   |
| Species Reactivity | Ms  |
| Specificity        | Natural and recombinant Mouse Arylsulfatase A |
| Target Name        | Arylsulfatase A                               |
| Other Names        | ASA, Cerebroside-sulfatase, As2               |
| Accession No.      | P50428  |

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.166ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Arsa

Note: This product is for in vitro research use only and is not intended for use in humans or animals.