

Mouse Anxa5/Annexin A5 ELISA Kit

Catalog No: #EK19255



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Anxa5/Annexin A5 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Annexin A5
Target Name	Annexin A5
Other Names	Anchorin CII,Annexin V,Annexin-5,Calphobindin I,CBP-I,Endonexin II,Lipocortin V,Placental anticoagulant protein 4,PP4,Placental anticoagulant protein I,PAP-I,Thromboplastin inhibitor,Vascular anticoagulant-alpha,VAC-alpha,Anx5
Accession No.	P48036

Application Details

Assay Type: SandwichDetect Range:1.56-100ng/mL Sensitivity:0.78ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Anxa5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.