

Pig B2M/Beta-2-microglobulin ELISA Kit

Catalog No: #EK19257



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pig B2M/Beta-2-microglobulin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Beta-2-microglobulin
Target Name	Beta-2-microglobulin
Other Names	Lactollin
Accession No.	Q07717

Application Details

Assay Type: Sandwich Detect Range: 3.12-200ng/mL Sensitivity: 1.62ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: B2M

Note: This product is for in vitro research use only and is not intended for use in humans or animals.