

Dog CTSK/Cathepsin K ELISA Kit

Catalog No: #EK19258



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Dog CTSK/Cathepsin K ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Cathepsin K
Target Name	Cathepsin K
Accession No.	Q3ZKN1

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.18ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:CTSK

Note: This product is for in vitro research use only and is not intended for use in humans or animals.