

# Rat Ctsk/Cathepsin K ELISA Kit

Catalog No: #EK19259



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Ctsk/Cathepsin K ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cathepsin K
Target Name	Cathepsin K
Accession No.	O35186

## Application Details

Assay Type: SandwichDetect Range:6.25-400pmol/L Sensitivity:3.125pmol/L Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:Ctsk

Note: This product is for in vitro research use only and is not intended for use in humans or animals.