

RBL2 (Phospho-Thr642) Antibody

Catalog No: #12696



Package Size: #12696-1 50ul #12696-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RBL2 (Phospho-Thr642) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	RBL2 (Phospho-Thr642) Antibody detects endogenous levels of RBL2 only when phosphorylated at Thr642
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RBL2 (Phospho-Thr642)
Target Name	RBL2
Modification	Phospho
Other Names	RBL2, p130, Retinoblastoma-like 2 (p130), RBR-2, Rb2, Retinoblastoma-like protein 2
Accession No.	Swiss-Prot#: Q08999NCBI Gene ID: 5934
Target Species	human
Calculated MW	128kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Images



Western Blot analysis of lysates of 1. HeLa cell and 2. HeLa cell Serum-free treated, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.